



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | G-129-C = G34120 |
| product description: | Human TnI Peptide 3 Polyclonal Antibody, Affinity Purified |
| lot number: | [lot specific] |
| immunogen: | Highly purified human cardiac troponin I |
| source: | Goat |
| form: | Liquid |
| purification method: | Antigen affinity purified (-RAYATEPHAKKKS-) |
| purity: | >90% (by SDS-PAGE) |
| buffer: | 10mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS) |
| preservative: | 0.05% (w/v) NaN ₃ |
| concentration: | [lot specific] mg/mL |
| expiration date: | [lot specific] |
| storage: | 2-8°C |
| country of origin: | Goat anti-serum from which this product is derived is of United States origin. |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |