



PRODUCT SPECIFICATION

CATALOG NUMBER: G-130-C

PRODUCT:	Goat anti-TNI peptide 4, affinity purified
IMMUNOGEN:	Highly purified human cardiac troponin I
FORM:	Liquid
PURITY:	>95% pure. Affinity purified with solid phase human cardiac troponin I peptide 4
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE:	2-8 degrees C
SPECIFICITY:	Cross reactivity with human skeletal Troponin I is less than 1% by ELISA.