



## PRODUCT SPECIFICATION

**CATALOG NUMBER: G-132-C**

PRODUCT:	Goat anti-TNI peptide 2, affinity purified
IMMUNOGEN:	Highly purified human cardiac Troponin I
FORM:	Liquid
PURITY:	Affinity purified
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE:	2-8 degrees C short term -20 degrees C or colder long term
SPECIFICITY:	Cross reactivity with human skeletal Troponin I is less than 1% by ELISA.