



PRODUCT SPECIFICATION

CATALOG NUMBER: G-149-C

PRODUCT:	Goat anti-GPBB peptide G, affinity pure
IMMUNOGEN:	Synthetic peptide (QIPPPNIPRD, residues 834-843) coupled to BSA
FORM:	Liquid
PURITY:	Affinity purified >95% by SDS-PAGE
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE:	2-8 degrees C