



PRODUCT SPECIFICATION

catalog number:	G-150-C = G34250015
product description:	TNI Peptide 6 Polyclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Highly purified human cardiac troponin I
source:	Goat
form:	Liquid
purification method:	Affinity purified
buffer:	10mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS)
preservative:	0.05% (w/v) NaN ₃
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C
country of origin:	Goat anti-serum from which this product is derived is of United States origin.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.