



PRODUCT SPECIFICATIONS

CATALOG NUMBER: G-202-C

PRODUCT:	Goat anti-mouse IgG Fc, affinity pure
FORM	Liquid
PURITY:	Affinity purified
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	This antibody reacts with heavy chains on mouse IgG, but not with light chains on mouse immunoglobulins. No reactivity is observed to non-immunoglobulin mouse serum proteins.