



PRODUCT SPECIFICATION

CATALOG NUMBER: G-301-C

PRODUCT:	Goat anti-rabbit IgG (H & L), affinity pure
FORM:	Liquid
PURITY:	Affinity purified
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE:	2-8 degrees C short term -20 degrees C or colder long term
SPECIFICITY:	This antibody reacts with heavy chains on rabbit IgG, and with light chains on all rabbit immunoglobulins. No reactivity is observed to non-immunoglobulin rabbit serum proteins.