



## PRODUCT SPECIFICATION

**CATALOG NUMBER: G-601-C**

PRODUCT:	Goat anti-rat IgG (H & L), affinity pure
FORM:	Liquid
PURITY:	Affinity purified.
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) Sodium Azide
STORAGE:	2-8°C
SPECIFICITY:	This antibody reacts with heavy chains on rat IgG and with light chains on all rat immunoglobulins by IEP. No reactivity is observed to non-immunoglobulin rat serum proteins by IEP.