



PRODUCT SPECIFICATION

catalog number:	G10800
product description:	NT-proBNP Polyclonal Antibody, Purified
lot number:	[lot specific]
source:	Goat
form:	Liquid
purification method:	Affinity purified
buffer:	10 mM sodium phosphate, 0.8% NaCl, pH 7.2; filtered in 0.2 µm solution
preservative:	0.05% NaN ₃
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for the 1-24aa sequence of NT-proBNP.
pairing:	This polyclonal can be used in ELISA as the detection antibody in conjunction with either monoclonal A10760 or A10770 (as the capture antibody).
manufacture date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.