



PRODUCT SPECIFICATION

CATALOG NUMBER: G12020

| | |
|----------------------|--|
| PRODUCT DESCRIPTION: | Goat Anti-LH Polyclonal, Affinity Purified |
| HOST ANIMAL: | Goat |
| FORM: | Liquid |
| BUFFER: | 0.05 M PBS, pH 7.4 |
| PRESERVATIVES: | 0.1% NaN ₃ |
| AFFINITY CONSTANT: | 3 x 10 ¹⁰ |
| STORAGE TEMPERATURE: | -20°C long term; 4-8°C short term |
| CROSSREACTIVITY: | < 0.01% FSH, < 0.01% TSH |
| PAIRING: | Pairs with LH Mab A20120 |
| HAZARD/BIOHAZARD: | This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |