



PRODUCT SPECIFICATION

catalog number:	G13010
product description:	HCG Polyclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human chorionic gonadotropin
source:	Rabbit
form:	Liquid
purification method:	Affinity purified
buffer:	0.05 M PBS, pH 7.4
preservative:	0.1% NaN ₃
concentration:	[lot specific] mg/mL (A _{280nm})
affinity constant:	9 x 10 ¹⁰
cross-reactivity:	▪ FSH [lot specific] ▪ LH [lot specific]
pairing:	This antibody can pair with hCG MAbs A35240, A35050 or A35230.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.