



## PRODUCT SPECIFICATION

catalog number:	<b>G13010</b>
product description:	HCG Polyclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human chorionic gonadotropin
source:	Rabbit
form:	Liquid
purification method:	Affinity purified
buffer:	0.05 M PBS, pH 7.4
preservative:	0.1% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL (A <sub>280nm</sub> )
affinity constant:	9 x 10 <sup>10</sup>
cross-reactivity:	▪ FSH [lot specific]                      ▪ LH [lot specific]
pairing:	This antibody can pair with hCG MAbs A35240, A35050 or A35230.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.