



PRODUCT SPECIFICATION

catalog number:	G14020126
product description:	HCG Alpha Polyclonal Antibody, Affinity Purified
lot number:	[lot specific]
immunogen:	Human alpha-HCG subunit
source:	Goat
form:	Liquid
purification method:	Affinity purified
buffer:	PBS, pH 7.4
preservative:	0.1% NaN ₃
concentration:	[lot specific] mg/mL
specificity:	This polyclonal antibody is alpha subunit specific and reacts with LH, hCG, FSH and Prolactin.
storage:	2-8°C; DO NOT FREEZE
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.