



PRODUCT SPECIFICATION

catalog number:	G19001
product description:	Human Azurocidin Polyclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant human azurocidin, NSO-derived (Ile27-Pro250)
source:	Goat
form:	Liquid
purification method:	Antigen affinity chromatography; filtered in 0.2 µm solution
buffer:	PBS, pH 7.4
preservative:	none
isotype:	IgG
concentration:	[lot specific] mg/mL
endotoxin:	<0.10 EU/µg by LAL
specificity:	This polyclonal antibody is specific for human azurocidin.
storage:	Upon receipt: -20°C (up to 12 months); After opening: 2-8°C (one month); -20°C (up to 6 months)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.