



PRODUCT SPECIFICATION

catalog number:	G44020
product description:	Human HBsAg Polyclonal Antibody, Affinity Purified
lot number:	[lot specific]
immunogen:	Purified human HBsAg
source:	Goat
form:	Liquid
purification method:	Affinity purified from whole goat serum
buffer:	PBS, pH 7.4
preservative:	0.02% NaN ₃
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.