



## PRODUCT SPECIFICATION

**CATALOG NUMBER: G47010**

PRODUCT DESCRIPTION:	Affinity Purified Rabbit anti- <u>Streptococcus A</u> , Polyclonal
HOST ANIMAL:	Rabbit
IMMUNOGEN:	Group A Streptococcus
PURIFICATION METHOD:	Affinity Chromatography
BUFFER:	PBS, pH 7.0
PRESERVATIVES:	None added
STORAGE:	4°-8C short term, -20°C long term
INTENDED USE:	For manufacturing use only