



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | G47020 |
| product description: | Streptococcus B Polyclonal Antibody, Affinity Purified |
| lot number: | [lot specific] |
| immunogen: | Streptococci, group B |
| source: | Rabbit |
| form: | Liquid |
| purification method: | Immunoaffinity chromatography |
| buffer: | 0.01M PBS, pH 7.2 |
| preservative: | 0.1% NaN ₃ |
| concentration: | [lot specific] mg/mL (A _{280nm}) |
| specificity: | This polyclonal antibody is type specific for the carbohydrate of group B Streptococcus. |
| cross-reactivity: | This antibody does not cross-react with other Strep groups. |
| storage: | 2-8°C short term (< 6 months); -20°C long term avoiding multiple freeze thaw cycles. |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |