



## PRODUCT SPECIFICATION

catalog number:	<b>G47090</b>
product description:	Streptococcus A Polyclonal Antibody, Affinity Purified
lot number:	[lot specific]
immunogen:	Streptococci, group A
source:	Goat
form:	Liquid
purification method:	Immunoaffinity chromatography
buffer:	0.01M phosphate buffer saline, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
concentration:	1.0 mg/mL
specificity:	Type specific carbohydrate of group A Streptococcus
cross-reactivity:	Does not cross react with other Strep groups
pairing:	This antibody can be used for both capture and detection in a pair.
notes/applications:	Conjugates well with colloidal gold.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.