



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>G47091</b>   |
| product description: | Streptococcus A Polyclonal Antibody, Affinity Purified  |
| lot number:          | [lot specific]  |
| immunogen:           | Streptococci, group A   |
| source:              | Rabbit  |
| form:                | Liquid  |
| purification method: | Immunoaffinity chromatography   |
| buffer:              | 0.01M phosphate buffer saline, pH 7.2   |
| preservative:        | 0.1% NaN <sub>3</sub>   |
| concentration:       | 1.0 mg/mL   |
| specificity:         | Type specific carbohydrate of group A Streptococcus   |
| cross-reactivity:    | Does not cross react with other Strep groups  |
| storage:             | 2-8°C   |
| hazard/biohazard:    | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |