



PRODUCT SPECIFICATION

CATALOG NUMBER: G49400

PRODUCT:	Cystatin C Polyclonal, Affinity Pure
HOST ANIMAL:	Goat
LOT NUMBER:	Reference
PURIFICATION:	Affinity purified using solid-phase Human Cystatin C
BUFFER:	10 mM sodium phosphate, 0.15 M Sodium Chloride, pH 7.2, 0.2 µm filtered
PRESERVATIVE:	0.05% NaN ₃
CONCENTRATION:	Reference (E ^{1%} at 280 nm = 13.0)
STORAGE CONDITIONS:	2-8°C
COUNTRY OF ORIGIN:	USA
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.