



PRODUCT SPECIFICATION

catalog number:	G49420
product description:	Human Cystatin C Polyclonal Antibody, Purified
lot number:	[lot specific]
source:	Rabbit
form:	Liquid
purification method:	Protein G chromatography
purity:	> 90%
buffer:	PBS, pH 7.4
preservative:	0.3% NaN ₃
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.