



PRODUCT SPECIFICATION

catalog number:	G56020496
product description:	Digoxigenin Polyclonal Antibody, Affinity Purified
lot number:	[lot specific]
immunogen:	Digoxigenin conjugated with KLH
source:	Goat
form:	Liquid
purification method:	Affinity purified
buffer:	10 mM PO ₄ , 150 mM NaCl, pH 7.4
preservative:	0.05% NaN ₃
concentration:	1.0 mg/mL
cross-reactivity:	This polyclonal antibody react 100% with digoxigenin and digoxin, it does not cross-react with carrier protein KLH or steroids
storage:	2-8°C; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.