



PRODUCT SPECIFICATION

catalog number:	G61401								
product description:	IGFBP-1 Polyclonal Antibody, Affinity Purified								
lot number:	[lot specific]								
immunogen:	Recombinant human IGFBP-1 (Ala26-Asn259), <i>E. coli</i> derived								
source:	Goat								
form:	Lyophilized								
purification method:	Antigen affinity purified								
buffer:	Lyophilized from a 0.2 µm filtered solution in PBS with trehalose								
preservative:	none								
isotype:	IgG								
quantity:	[lot specific] mg								
endotoxin:	<0.10 EU/µg by LAL								
specificity:	This polyclonal antibody is specific for human IGFBP-1.								
cross-reactivity:	<table><tr><td>▪ rh IGFBP-2</td><td>< 2.0%</td><td>▪ rh IGFBP-4</td><td>< 2.0%</td></tr><tr><td>▪ rh IGFBP-3</td><td>< 2.0%</td><td>▪ rh IGFBP-5</td><td>< 2.0%</td></tr></table>	▪ rh IGFBP-2	< 2.0%	▪ rh IGFBP-4	< 2.0%	▪ rh IGFBP-3	< 2.0%	▪ rh IGFBP-5	< 2.0%
▪ rh IGFBP-2	< 2.0%	▪ rh IGFBP-4	< 2.0%						
▪ rh IGFBP-3	< 2.0%	▪ rh IGFBP-5	< 2.0%						
expiration date:	[lot specific]								
reconstitution:	Reconstitute to 0.2mg/mL with sterile PBS.								
storage:	Lyophilized: shipped at 2-8°C; upon receipt store -20°C (up to 12 months) After reconstitution: 2-8°C (one month); -20°C (up to 6 months) avoiding multiple freeze thaw cycles.								
hazard/biohazard:	There is no known biohazard associated with this product.								