



PRODUCT SPECIFICATION

CATALOG NUMBER: G62020

PRODUCT DESCRIPTION:	Goat Anti-Human CEA, Affinity Purified
HOST ANIMAL:	Goat
FORM:	Liquid
PURIFICATION:	CEA affinity column chromatography purified immunoglobulin fractions
BUFFER:	0.015 M KPO ₄ , 0.85% NaCl, pH 7.4
PRESERVATIVES:	0.1% NaN ₃
STORAGE TEMPERATURE:	2-8°C
CROSS-REACTIVITY:	No cross-reactivity with human albumin , hCG, AFP, Ferritin, PSA, and PAP
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.