



PRODUCT SPECIFICATION

catalog number:	G83460 = D-401-C-ABS
product description:	Goat IgG (H&L) Antibody, Purified
lot number:	[lot specific]
source:	Donkey
form:	Liquid
purification method:	Affinity purification
buffer:	10mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS)
preservative:	0.05% NaN ₃
concentration:	[lot specific] mg/mL
specificity:	This antibody reacts with goat IgG.
cross-reactivity:	This antibody does not react to human IgG, rabbit IgG, mouse IgG, or rat IgG.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.