



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | G83460 = D-401-C-ABS |
| product description: | Goat IgG (H&L) Antibody, Purified |
| lot number: | [lot specific] |
| source: | Donkey |
| form: | Liquid |
| purification method: | Affinity purification |
| buffer: | 10mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS) |
| preservative: | 0.05% NaN ₃ |
| concentration: | [lot specific] mg/mL |
| specificity: | This antibody reacts with goat IgG. |
| cross-reactivity: | This antibody does not react to human IgG, rabbit IgG, mouse IgG, or rat IgG. |
| expiration date: | [lot specific] |
| storage: | 2-8°C |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |