



PRODUCT SPECIFICATION

catalog number:	G83720
product description:	Fibrinogen Polyclonal Antibody, Affinity Purified
lot number:	[lot specific]
immunogen:	Human fibrinogen
source:	Goat
form:	Liquid, 0.2 μ m filtered
purification method:	Antigen affinity chromatography
purity:	\geq 95% by SDS-PAGE
buffer:	10mM sodium phosphate, 0.15M NaCl, pH 7.2
preservative:	0.05% (w/v) NaN ₃
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 13.0)
specificity:	This antibody reacts with human fibrinogen.
expiration date:	[lot specific]
storage:	2-8°C
country of origin:	U.S.A.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.