



## PRODUCT SPECIFICATION

catalog number:	<b>G91000</b>
product description:	PCT (N-terminal) Polyclonal Antibody
lot number:	[lot specific]
source:	Goat
form:	Liquid
purification method:	Antigen affinity purified, 0.2 $\mu$ filtered
purity:	>85% by SDS-PAGE
buffer:	10mM Na <sub>3</sub> PO <sub>4</sub> , 0.8% NaCl, pH 7.2
preservative:	0.05% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL
manufacture date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.