



## PRODUCT SPECIFICATION

catalog number:	<b>G91010</b>
product description:	Procalcitonin (PCT) N-terminal Polyclonal Antibody, Affinity Purified
lot number:	[lot specific]
source:	Goat
form:	Liquid
purification method:	Antigen affinity purified, 0.2 $\mu$ filtered
purity:	>85% by SDS-PAGE
buffer:	10mM sodium phosphate, 0.8% NaCl, pH 7.2
preservative:	0.05% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL
pairing:	This antibody can be used in ELISA as the detection antibody in conjunction with A91082 (as the capture antibody).
manufacture date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.