



PRODUCT SPECIFICATION

catalog number:	G91040
product description:	Human Procalcitonin Polyclonal Antibody, Affinity Purified
lot number:	[lot specific]
immunogen:	Recombinant human Procalcitonin (Ala26-Gly117), E. coli-derived
source:	Sheep
form:	Lyophilized
purification method:	Antigen affinity purified; 0.2µm filtered
buffer:	PBS with Trehalose, pH 7.4
preservative:	None
isotype:	IgG
concentration:	[lot specific] mg/mL
specificity:	This polyclonal antibody is specific for human Procalcitonin. .
reconstitution:	Reconstitute at 0.2 mg/mL in sterile PBS
storage:	Upon receipt store -20°C (up to 12 months); After opening: 2-8°C (one month); -20°C (up to 6 months)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.