



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: G91200**

PRODUCT DESCRIPTION:	Goat Anti-GC Globulin (Vitamin D Binding Protein), Affinity Purified
FORM:	Liquid, 0.2 µm filtered
PURIFICATION:	Affinity purified using solid phase human GC Globulin
BUFFER:	10 mM sodium phosphate, 0.15 M sodium chloride, pH 7.2
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE TEMPERATURE:	2-8° C
COUNTRY OF ORIGIN:	USA
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution.