



PRODUCT SPECIFICATION

CATALOG NUMBER: J02080

PRODUCT DESCRIPTION:	Dengue Virus NS1, C terminal antigen, recombinant
SOURCE:	<u>E. coli</u> , recombinant
PURITY:	>90% pure by SDS-PAGE
BUFFER:	8M urea, 20mM Tris-HCl, 60 mM NaCl, pH 8.0
FUSION:	None
STORAGE TEMPERATURE:	2-8°C short term -20°C long term
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.