



PRODUCT SPECIFICATION

CATALOG NUMBER: J02090

PRODUCT DESCRIPTION:	Dengue Virus NS1 Antigen, N terminus epitope, recombinant
SOURCE:	<i>E. coli</i>
PURITY:	>90% pure by SDS-PAGE
FUSION:	Glutathione transferase
BUFFER:	100mM Tris-HCl, 60 mM NaCl, 5% Glycerol, pH 8.0
STORAGE:	2-8°C short term; -20°C long term
APPLICATION:	ELISA, WB
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.