



PRODUCT SPECIFICATION

catalog number:	J030005
product description:	SARS-CoV-2 Spike S1 Recombinant Antigen, aa 14-681
product description:	Protein contains extra amino acids (AS) in N-terminus and His-6 tag at C-terminus in the protein.
lot number:	[lot specific]
source:	CHO cells; recombinant
form:	Liquid
purity:	>95% by SDS-PAGE
quality control:	Coomassie strained SDS-PAGE, analytical SEC
buffer:	PBS, pH 7.4
preservative:	None
molecular weight:	76.42 kDa
concentration:	[lot specific] mg/mL
specificity:	Immunoreactive with sera of COVID-19 infected individuals. Binds with ACE2.
manufacture date:	[lot specific]
storage:	-20°C
hazard/biohazard:	There is no known hazard or biohazard associated with this product.