



## PRODUCT SPECIFICATION

catalog number:	<b>J03007</b>
product description:	SARS-CoV-2 S1 Nucleoprotein Recombinant Antigen, Full Length
lot number:	[lot specific]
source:	CHO cells; recombinant
form:	Liquid
purity:	>90% by SDS-PAGE
quality control:	Coomassie strained SDS-PAGE, analytical SEC
buffer:	PBS, pH 7.4
preservative:	None
molecular weight:	46.9 kDa
concentration:	[lot specific] mg/mL
specificity:	Immunoreactive with sera of COVID-19 infected individuals.
manufacture date:	[lot specific]
storage:	-20°C
hazard/biohazard:	There is no known hazard or biohazard associated with this product.