



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J03010**

PRODUCT DESCRIPTION:	SARS-Associated Coronavirus Spike mosaic S(N), recombinant antigen
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	Sepharose derived purification
PURITY:	>95% pure by 10% PAGE (Coomassie staining)
BUFFER:	25 mM Tris HCl, 0.4% sarcosyl, 0.25% Triton-100, 50% glycerol
STORAGE TEMPERATURE:	Shipped at ambient temperature. Upon arrival, store at -20°C.
STABILITY:	Six months at 2-8°C, five years at -20°C.
SPECIFICITY:	Immunoreactive with sera of SARS-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.