



PRODUCT SPECIFICATION

catalog number:	J03014
product description:	SARS-CoV-2 Delta Variant Spike Protein, Recombinant, with C-terminal 6-His Tag Delta variant B.1.617.2 T19R, G142D, F157del, R158del, L452R, T478K, D614G, P681R, D950N
N-terminal sequence	Val16
lot number:	[lot specific]
source:	HEK293
form:	Lyophilized from a 0.2um filtered solution in PBS with Trehalose
purity:	>95% by SDS-PAGE
buffer:	PBS
preservative:	None
molecular weight:	138 kDa
concentration:	[lot specific] mg/ml
endotoxin:	<0.10 EU/1ug by LAL method
reconstitution:	Reconstitute at 500ug/ml in PBS
storage:	Upon receipt; -20°C to -70°C (up to 12 months); 2-8°C (one month after opening); -20°C to -70°C (up to 3 months after opening)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.