



PRODUCT SPECIFICATION

catalog number:	J03019
product description:	SARS-CoV-2 Delta Variant Spike s1 Recombinant Antigen.
product information:	Protein contains amino acids 14-681, mutations T19R, G142D, E156G, del_F157-R158 L452R, T478K, D614G, P681R, two extra amino acids (AS) in N-terminus and His-6 tag at C-terminus and GSG linker between protein and tag.
lot number:	[lot specific]
source:	CHO cells; recombinant
form:	Liquid
purity:	>95% by SDS-PAGE
quality control:	SDS-PAGE, analytical SEC, Octet binding
buffer:	PBS, pH 7.4
preservative:	None
concentration:	[lot specific] mg/mL
manufacture date:	[lot specific]
expiration date:	One year from date of receipt
storage:	-70°C
hazard/biohazard:	There is no known hazard or biohazard associated with this product.