



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J03030**

PRODUCT DESCRIPTION:	Recombinant SARS-Associated Coronavirus Spike mosaic (S/C)
PURIFICATION METHOD:	Purified by chromatography.
PURITY:	>95% determined by 10% PAGE (coomassie staining)
BUFFER:	25 mM Tris-HCl, 0.4% sarcosyl, 0.25% Triton-100 and 50% glycerol
STORAGE TEMPERATURE:	-20°C
STABILITY:	Five years frozen, 6 months at +4°C