



PRODUCT SPECIFICATION

catalog number:	J03065
product description:	SARS-CoV-2 Spike Glycoprotein (S) (full sequence), Recombinant Antigen, mammalian expression
lot number:	[lot specific]
source:	Hek293 cells
form:	Liquid
purity:	>90%
buffer:	Tris buffer with additives and salts, pH 8.0
concentration:	[lot specific]mg/ml
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.