



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J03070**

PRODUCT DESCRIPTION:	Recombinant SARS Associated Coronavirus M protein.
SOURCE:	<u>E. coli</u> expressed recombinant with GST fusion tag.
PURIFICATION METHOD:	Purified by affinity chromatography.
PURITY:	>95% determined by SDS-PAGE
BUFFER:	50 mM Tris-HCl, 60 mM NaCl, with 50% Glycerol.
STORAGE TEMPERATURE:	-20°C
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.