



## PRODUCT SPECIFICATION

catalog number:	J03073
product description:	OC43 Coronavirus Nucleoprotein (NP)(full sequence) Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E.Coli</i>
form:	Liquid
purity:	>90%
buffer:	Tris buffer with additives and salts, pH 8.0
concentration:	[lot specific] mg/ml
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.