



PRODUCT SPECIFICATION

CATALOG NUMBER: J03080

PRODUCT DESCRIPTION:	Recombinant SARS Associated Coronavirus E protein.
SOURCE:	<u>E. coli</u> expressed recombinant with GST fusion tag.
PURIFICATION METHOD:	Purified by affinity chromatography.
PURITY:	>95% determined by SDS-PAGE
BUFFER:	50 mM Tris-HCl, 60 mM NaCl, with 50% Glycerol.
STORAGE TEMPERATURE:	Shipped at 2-8°C Store at -20°C
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.