



PRODUCT SPECIFICATION

catalog number:	J03083
product description:	SARS-Associated Coronavirus N (Nucleocapsid) Recombinant Antigen
lot number:	[lot specific]
source:	E. coli
form:	Liquid
purity:	90% by SDS-PAGE
buffer:	PBS
concentration:	1 mg/ml
specificity:	aa 1 to 49 of the N protein immunodominant region
applications:	ELISA, WB
storage:	2-8°C short term (<3 months); -80°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.