



PRODUCT SPECIFICATION

catalog number:	J03084
product description:	SARS-Associated Coronavirus Nucleocapsid Recombinant Antigen
lot number:	[lot specific]
source:	E. coli
form:	Liquid
purity:	90% by SDS-PAGE
buffer:	1.5 M Urea; 20mM Tris-HCL, pH 7.2 / 50% Glycerol
concentration:	1 mg/ml
specificity:	aa 340-390 immunodominant region
applications:	ELISA, WB
storage:	2-8°C short term (<3 months); -80°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.