



PRODUCT SPECIFICATION

catalog number:	J03088
product description:	SARS-Associated Coronavirus Spike Mosaic S(C) Recombinant Antigen
lot number:	[lot specific]
source:	E. coli
form:	Liquid
purity:	90% by SDS-PAGE
buffer:	1.5M Urea; 20mM Tris-HCl, pH 7.2 / 50% Glycerol
concentration:	1 mg/ml
specificity:	Spike protein immunodominant regions
applications:	ELISA, WB
storage:	2-8° short term (<3 months); -80°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.