



PRODUCT SPECIFICATION

catalog number:	J03089
product description:	Coronavirus-19 Spike-E-M Mosaic Recombinant Antigen, immunodominant regions
lot number:	[lot specific]
source:	E. coli
fusions:	Six histidines at C-terminus
form:	Liquid
purity:	>90% by SDS-PAGE
buffer:	1 x PBS
concentration:	1 mg/ml
specificity:	aa 20-210 Spike protein immunodominant regions; aa 51-75 E immunodominant regions; aa 100-222 M immunodominant regions
cross reactivity:	Reactivity with SARS infected individuals not tested.
applications:	ELISA, WB
storage:	2-8°C short term (<3 months); -80°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.