



## PRODUCT SPECIFICATION

catalog number:	J03093
product description:	Coronavirus-19 Spike Receptor Binding Domain (aa 300-600) Recombinant Antigen
lot number:	[lot specific]
source:	<i>E. coli</i>
fusions:	Six histidines at C-terminus
form:	Liquid
purity:	>90% by SDS-PAGE
buffer:	PBS
concentration:	1 mg/ml
specificity:	aa 300 to 600 of the Spike immunodominant region. Reactivity with SARS-CoV-2 infected individuals not tested.
applications:	ELISA, WB
storage:	2-8°C short term (<3 months); -80°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.