



PRODUCT SPECIFICATION

catalog number:	J10145566
product description:	Human IL-1 R4 (ST2) Antigen, Recombinant
lot number:	[lot specific]
source:	NSO derived, Lys19-Phe328 with C-terminal 6-His tag; N-term is Lys19
form:	Liquid; 0.2µm filtered
purity:	>97% by SDS-PAGE
buffer:	PBS, pH 7.4
molecular weight:	~100 kDa
concentration:	[lot specific] mg/mL
endotoxin level:	<0.10 EU/1µg by LAL
manufacture date:	[lot specific]
storage:	Upon receipt store -20°C (up to 12 months); After opening: 2-8°C (one month); -20°C (up to 6 months)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.