



PRODUCT SPECIFICATION

CATALOG NUMBER: J10270

PRODUCT DESCRIPTION:	Interleukin 18 (IL-18) Recombinant Antigen, High Pure
SOURCE:	Bacterial
FORM:	Liquid
PURITY RATIO:	>95% by SDS-PAGE Single band @ 18.2 kD
BUFFER:	10 mM phosphate buffer saline, 20% glycerol
IL-18 CONTENT:	>1 mg/ml (by ELISA, MBL Japan)
BIO-BURDEN:	<10 CFU/ml
STORAGE TEMPERATURE:	-70°C. Avoid multiple freeze thaw cycles.
HAZARD/BIOHAZARD:	All donor units have been tested and found negative for HBsAg, anti-HCV and anti-HIV I & II antibodies by FDA approved methods; found negative for viral HBV, HCV and HIV I by PCR.